
Contents

<i>Preface</i>	<i>v</i>
<i>Contributors</i>	<i>ix</i>
1 New Insights into GPCR Function: Implications for HTS	1
<i>Richard M. Eglén and Terry Reisine</i>	
2 Screening Technologies for G Protein-Coupled Receptors: From HTS to uHTS	15
<i>Maite de los Frailes and Emilio Díez</i>	
3 GPCR Signaling: Understanding the Pathway to Successful Drug Discovery	39
<i>Christine Williams and Stephen J. Hill</i>	
4 An Overview on GPCRs and Drug Discovery: Structure-Based Drug Design and Structural Biology on GPCRs	51
<i>Kenneth Lundström</i>	
5 Understanding the Ligand–Receptor–G Protein Ternary Complex for GPCR Drug Discovery	67
<i>Venkata R.P. Ratnala and Brian Kobilka</i>	
6 Assay Data Quality Assessment	79
<i>Hanspeter Gubler</i>	
7 Homology Modeling of GPCRs	97
<i>John Simms, Nathan E. Hall, Polo H.C. Lam, Laurence J. Miller, Arthur Christopoulos, Ruben Abagyan, and Patrick M. Sexton</i>	
8 GPCR Expression Using Baculovirus-Infected Sf9 Cells	115
<i>Amanda L. Aloia, Richard V. Glatz, Edward J. McMurchie, and Wayne R. Leifert</i>	
9 Radioligand Binding Assays: Application of [¹²⁵ I]Angiotensin II Receptor Binding	131
<i>Wayne R. Leifert, Olgatina Bucco, Mahinda Y. Abeywardena, and Glen S. Patten</i>	
10 [³⁵ S]GTPγS Binding in G Protein-Coupled Receptor Assays	143
<i>Tamara Cooper, Edward J. McMurchie, and Wayne R. Leifert</i>	
11 A Time-Resolved Fluorescent Lanthanide (Eu)-GTP Binding Assay for Chemokine Receptors as Targets in Drug Discovery	153
<i>Jean Labrecque, Rebecca S.Y. Wong, and Simon P. Fricker</i>	
12 Use of the DiscoverX Hithunter cAMPII Assay for Direct Measurement of cAMP in Gs and Gi GPCRs	171
<i>Joe Bradley and David McLoughlin</i>	
13 Use of Aequorin for G protein-Coupled Receptor Hit Identification and Compound Profiling	181
<i>Stephen J. Brough and Parita Shah</i>	

14	BacMam: Versatile Gene Delivery Technology for GPCR Assays	199
	<i>Elizabeth A. Davenport, Parvathi Nuthulaganti, and Robert S. Ames</i>	
15	Yeast Assays for G Protein-Coupled Receptors	213
	<i>Simon J. Dowell and Andrew J. Brown</i>	
16	GPCR Microspot Assays on Solid Substrates	231
	<i>Ye Fang and Joydeep Lahiri</i>	
17	Resonant Waveguide Grating Biosensor for Whole-Cell GPCR Assays	239
	<i>Ye Fang, Ann M. Ferrie, and Elizabeth Tran</i>	
18	FRET-Based Measurement of GPCR Conformational Changes	253
	<i>Sébastien Granier, Samuel Kim, Juan José Fung, Michael P. Bokoch, and Charles Parnot</i>	
19	FLIPR [®] Assays of Intracellular Calcium in GPCR Drug Discovery	269
	<i>Kasper B. Hansen and Hans Bräuner-Osborne</i>	
20	Use of Fluorescence Indicators in Receptor Ligands	279
	<i>Kaleeckal G. Harikumar and Laurence J. Miller</i>	
21	Application of Fluorescence Resonance Energy Transfer Techniques to Establish Ligand–Receptor Orientation	293
	<i>Kaleeckal G. Harikumar and Laurence J. Miller</i>	
22	Detection of GPCR/ β -Arrestin Interactions in Live Cells Using Bioluminescence Resonance Energy Transfer Technology	305
	<i>Martina Kocan and Kevin D.G. Pflieger</i>	
23	Using Reporter Gene Technologies to Detect Changes in cAMP as a Result of GPCR Activation.	319
	<i>Daniel J. Rodrigues and David McLoughlin</i>	
24	A Quantum Dot-Labeled Ligand–Receptor Binding Assay for G Protein-Coupled Receptors Contained in Minimally Purified Membrane Nanopatches	329
	<i>Jody L. Swift, Melanie C. Burger, and David T. Cramb</i>	
25	<i>Xenopus</i> Oocyte Electrophysiology in GPCR Drug Discovery	343
	<i>Kasper B. Hansen and Hans Bräuner-Osborne</i>	
26	Immunoprecipitation and Phosphorylation of G Protein-Coupled Receptors	359
	<i>Walter G. Thomas</i>	
	<i>Index</i>	373